

## Thermo Pal

Medical **Non-contact** Thermometer

**Medical Non-contact thermomter. Medizinsische kontaktlose thermometer.Termómetro clínico sin contacto.Thermomètre médical sans contact .Non-contact termometru medical**

- Display
- Display
- Pantalla
- Écran
- Afişaj

- Measuring sensor
- Messsensor
- Sensor de medición
- Capteur de mesure
- Senzor de măsurare

- Battery compartment
- Batteriefach
- Compartimiento de la batería
- Compartiment de la batterie
- Compartiment pentru baterii

EC REF

Name: Niish Technologies Ltd.  
Address: Office 8, Marcus House, Park Hall Business Village,Longton, Stoke-on-Trent, ST2 5XA, United Kingdom.

REF

ETH-301

digiO2 International Co., Ltd.

website: <http://www.digio2.com/>

email: [service@digio2.com](mailto:service@digio2.com)


3F., No.582, Guohua Rd., Miaoli City, Miaoli County 360, Taiwan

© 201X digiO2 International Co., Ltd. All rights reserved.

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EN

Thank you for purchasing Thermo Pal medical **non-contact** thermometer. Please read this instruction manual carefully prior to first time use.

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### FEATURES

- Measures without skin contact (in 5cm from the forehead)
- Easy to read : 2 inches large screen
- Quick Check \*: 3-color LCD backlight for temperature indication
- Long-Term Tracking \*: 10-memory recall with time and date display
- Communicationable \*:Bluetooth data transmission
- \* Depends on models

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### IMPORTANT INFORMATION

- Clean the device with a dry, soft cloth or cloth moistened in alcohol. Never use thinner, benzine or cleaner with abrasives.
- This device is not water resistant. Protect it from liquid spills.
- Keep this device out of the reach of children. A child may swallow the battery. If a child swallows them, seek medical treatment immediately.
- Do not self-diagnose your condition using the measured result. Consult your doctor if your temperature is higher than the normal temperature or you feel unwell.
- Use the thermometer in a room temperature between 10°C and 40°C.
- Always use the thermometer with its probe tip clean.
- To get an accurate reading, please push hair aside and wipe perspiration from the forehead when measuring.

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CE  
0120

= Type BF

Specifications are subject to change without prior notice  
All brand names and product names are trademarks or registered trademarks of their respective owners.

Always consult your doctor. Do not attempt self-diagnosis or self-treatment based on the measurement results and analysis. Self-diagnosis or self-treatment may lead to deterioration in your condition.

EN

OPERATION

- Measure body temperature: Press "MENU" or "MEASURE" to turn on the device. Press "MENU" to choose modes and press "MEASURE" to take measurement.

	LCD Backlights*	Symbol
36.0°C~37.8°C	Green	
34°C~35.9°C	Yellow	
37.9°C~43°C	Red	

\* Depends on models.

- Check records :  
Hold "MENU" in 2 seconds to check records in turn. The latest record is shown at first and press "MENU" to check other records.
- Change unit \*  
Hold "MEASURE" button until batteries are inserted to change °C and °F.
- Bluetooth transmission \*  
The data will be transferred automatically after the device displays the results.  
\* Depends on models.

TROUBLE SHOOTING

- The screen displays "LO":  
The measured temperature is less than 34 °C (93.2 °F)
- The screen displays "HI":  
The measured temperature is greater than 43 °C (109.4 °F)
- The screen displays the battery indicator " " and sounds 1 long and 2 short beeps:  
Please insert a set of new batteries.
- Blank screen when the power is turned on
  1. Please check battery position.
  2. Please replace new batteries.

Marks and Abbreviations

- CE -European Conformity. This device conforms to the standards for products sold in the European Economic Area.
- RoHS - This device conforms to the European Union Restriction of Hazardous Substances Directive.

- Caution: Do not allow to become damp or immerse in water.
- In accordance with Waste Electrical and Electronic Equipment Directive (WEEE) this unit contains electronic parts and must be disposed of in accordance with local laws regarding electronic waste products.

- Manufacturer
- Catalogue number
- Serial number
- Batch code
- Operator's manual
- Representative in the European Community
- Type BF Equipment
- IP22 Ingress protection class

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digiO2 International Co., Ltd.

CE 0120 RoHS Type BF  
Specifications are subject to change without prior notice.  
All brand names and product names are trademarks or registered trademarks of their respective owners.

Always consult your doctor. Do not attempt self-diagnosis or treatment based on the measurement results and analysis. Self-diagnosis or self-treatment may lead to deterioration in your condition.

## Federal Communication Commission interference statement

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

**This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:**

- (1) this device may not cause harmful interference and**
- (2) this device must accept any interference received, including interference that may cause undesired operation**

## RF Exposure Warning

The equipment complies with RF exposure limits set forth for an uncontrolled environment. The antenna(s) used for this transmitter must not be co-located or operating in conjunction with any other antenna or transmitter.

You are cautioned that changes or modifications not expressly approved by the party responsible for compliance could void your authority to operate the equipment.

This equipment should be installed and operated with a minimum distance of 20 centimeters between the radiator and your body.